

A World of Difference

The Clinicians Programme

Aims and Objectives

Parts 1, 2 & 3

By the end of the programme the participants should better understand:

- the underlying principles of prosthetics
- how, by involving the patient, the clinical outcome can be improved
- how to assess the front tooth position
- how to improve their understanding of impression taking including primary, the French Impression, the upper working impression and the lower primary impression
- how to better optimise the border seal
- when and how to use magnetic attachments
- how to improve their knowledge of record rims, measuring the vertical dimension of the face and checking speech
- the purpose of the gothic arch technique and how to apply it
- how to check the anterior set up
- how to take the retruded (centric) jaw relationship
- how to take a face bow record
- how to choose teeth and gingival colouration
- how to use the follow-up visits to improve patient satisfaction

Approved for 5 hours
verifiable CPD*



The Technicians Programme

Aims and Objectives

Parts 1 & 2

By the end of the programme the participants should better understand:

- how to use patient photographs or casts taken from a relative to choose the correct width and shape of tooth
- the underlying principles of prosthetics in the dental laboratory
- how to set-up the anterior teeth
- how to set-up the posterior teeth
- the different techniques for waxing in natural features including S curves and the rugae and creating more natural surface effects
- the flasking of the denture and in-flask correction
- the principles of natural gingival colouring and how to reproduce it in the denture
- how to de-flask the denture so that the natural effects that have been created are not lost in doing so

Approved for 2.5 hours
verifiable CPD*



A World of Difference

A World of Difference

Schottlander have produced three films to inspire and inform the Clinician, Technician and Patient. These films show the work and teaching of Dr John Besford and Mrs Ruth Bourke who also advised on the development of the Enigma Denture System. These programmes are intended to give you a better understanding of the underlying principles of prosthetics, and to inspire you to create the most natural looking dentures.

The Clinicians Programme Parts 1, 2 & 3

In three parts totalling 80 minutes, gives tips and hints on a wide range of clinical procedures.

Approved for 5 hours verifiable CPD*



VHS

containing 2 videos

Code 0907-V

CD ROM*

containing 3 CD's, Approved for 5 hours verifiable CPD

Code 0907-C

The Technicians Programme Parts 1 & 2

In two parts totalling 40 minutes, explains how laboratory techniques can be adapted to today's cosmetic world.

Approved for 2.5 hours verifiable CPD*



VHS

containing 2 videos

Code 0908-V

CD ROM*

containing 2 CD's, Approved for 2.5 hours verifiable CPD

Code 0908-C

A World of Difference Box Set

This DVD contains the Clinicians Programmes, Parts 1, 2 & 3 and the Technicians Programmes, Parts 1 & 2.



DVD

This DVD contains the Clinicians Programmes, Parts 1, 2 & 3 and the Technicians Programmes, Parts 1 & 2

Code 0909

The Patients Programme

This programme gives patients a better understanding of the stages involved in making dentures including many that they would never normally see.

Through other patients sharing their own experiences on screen they are led to see the vital role they themselves must play in order for the Clinician to provide them with dentures that function better and look more natural.

The patients programme is 15 minutes long and is intended either for the waiting room or to be given to patients to take home.



SINGLE

VHS
CD Rom
DVD

Code 0906-V
Code 0906-C
Code 0906-D

MULTI-PACK

Pack of 5

VHS
CD Rom
DVD

Code 0906-VMP
Code 0906-CMP
Code 0906-DMP